

E 926
Module
for the optimum loading of long MIS instruments.



- Holds up to 13 hollow instruments

EAN: 4002514504353 / Material number: 07297050 / Old Material Number: 69592601D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert/module for MIS instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Capacity	
Lumened instruments [number]	13
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	40
External dimensions, net width in mm	568
External dimensions, net depth in mm	510
External dimensions, gross height in mm	280
External dimensions, gross width in mm	580
External dimensions, gross depth in mm	630
Net weight in kg	2,40
Gross weight in kg	2,60
Companion pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	5
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	2
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
E 456 spring stay for MIS instruments [number]	2
E 464 holder for injector nozzle E 454 [number]	1