

E 429/1

Mobile unit

for optimum cleaning of micro instruments (eye OP), for 2-4 OP sets.



- With integrated magnetic spray arm
- Level 1: e.g. 4 x insert E441/1 or 2 x insert E142 per level
- Level 2 left: 20 connections for hollow instruments
- Level 2 right: 16 Luer lock adapters, male, horizontal Hose
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002514556253 / Material number: 07333120 / Old Material Number: 69542902D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
Product category	
Mobile unit for ophthalmology instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of ophthalmology instruments	•
Capacity	
Theatre sets for micro instruments [number]	2-4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	520
External dimensions, net width in mm	540
External dimensions, net depth in mm	540
External dimensions, gross height in mm	530
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	2,40
Gross weight in kg	3,30
Companion pack	
E 476 holder [number]	1
E 477 stopper [number]	1
E 790 Luer Lock connection element female to female [number]	8
E 791 Luer Lock connection element female to male [number]	8