

E 931

Insert

for optimum loading of up to 20 theatre shoes up to size 48.



- For use in upper and lower levels of mobile unit E 975/2
- Equipped with 20 holders
- Total capacity per cycle: 40 theatre shoes (20 pairs)

EAN: 4002514598444 / Material number: 07439490 /
Old Material Number: 69593101D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert for theatre shoes	•
Area of application	
Reprocessing of theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Theatre shoes [number]	20
Standard electrical connection	
Number of phases	0
Dimensions and weight	
External dimensions, net height in mm	254
External dimensions, net width in mm	586
External dimensions, net depth in mm	752
External dimensions, gross height in mm	275
External dimensions, gross width in mm	605
External dimensions, gross depth in mm	765
Net weight in kg	4,35
Gross weight in kg	4,71