

E 944/2

Module

for the optimum loading of narrow necked glassware 500 - 1,000 ml.



- 15 nozzles (E 352) 6 x 220 mm with clips (E 354)

EAN: 4002514605050 / Material number: 07459410 / Old Material Number: 69594403D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	15
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	250
External dimensions, net width in mm	558
External dimensions, net depth in mm	352
External dimensions, gross height in mm	180
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	570
Net weight in kg	3,92
Gross weight in kg	4,42