

E 555/2

Mobile unit

for cleaning instruments on 2 levels.



- Loading volume from below: H 210, W 487, D 592 mm
- Loading volume from above: H 219, W 495, D 548 mm
- Magnetic strip for mobile unit recognition (without magnets ML)
- Wide range of fitments available

EAN: 4002514606248 / Material number: 07460290 / Old Material Number: 69555503D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
<b>Product category</b>	
Mobile unit for inserts	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of containers	•
<b>Capacity</b>	
DIN mesh trays [number]	4
SPRI mesh trays [number]	2
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	298
External dimensions, net width in mm	530
External dimensions, net depth in mm	645
External dimensions, gross height in mm	650
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	640
Net weight in kg	10,00
Gross weight in kg	12,00
<b>Companion pack</b>	
Magnetic strip for automatic mobile unit recognition (without magnets) [number]	1