

DOS S 20  
 Bellows dispenser pump  
 For dispensing liquid acid media.



- For dispensing rinsing agent, neutralising agent
- For retro-fitting by Miele Service

EAN: 4002514641904 / Material number: 07514600 / Old Material Number: 69226004D

Product classification	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8536	•
Product category	
Dispenser system	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	245
External dimensions, gross width in mm	258
External dimensions, gross depth in mm	270
Net weight in kg	2,45
Gross weight in kg	2,91