

RK LWMM  
 Conductivity meter  
 for retrofitting on PG 8536.



EAN: 4002514641928 / Material number: 07514620 / Old Material Number: 69224015D

Product classification	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8536	•
Product category	
Other accessories	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	280
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	420
Net weight in kg	2,02
Gross weight in kg	2,97