


E 491

Insert

for the optimum loading of rotating instruments.



- Suitable for 18 turbine shafts
- Suitable for 30 contra angle shafts
- Also suitable for rotating instruments, such as used in podiatry

EAN: 4002514697949 / Material number: 07600210 / Old Material Number: 68549101D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Processing of dental instruments	•
<b>Capacity</b>	
Sheaths [number]	48
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	53
External dimensions, net width in mm	70
External dimensions, net depth in mm	120
External dimensions, gross height in mm	100
External dimensions, gross width in mm	265
External dimensions, gross depth in mm	180
Net weight in kg	0,85
Gross weight in kg	1,00