

E 497

Nozzle

for optimum loading of double and single tube breathing masks. Height 195 mm.

EAN: 4002514811253 / Material number: 07759470 / Old Material Number: 69749701D



- H 195 mm, with pipe 109 mm at the side

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	195
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Net weight in kg	0,06
Gross weight in kg	0,07