

E 499

Rinse fittings

for the optimum loading of suction, rinse and Veress cannulae.



- Different connection types for suction & rinsing cannulae available
- Compatible with transfer trolleys with injector nozzles

EAN: 4002514811482 / Material number: 07759500 / Old Material Number: 69549901D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Irrigation sleeve	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
Reprocessing of ENT instruments	•
Reprocessing of lumened instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	30
External dimensions, gross width in mm	110
External dimensions, gross depth in mm	170
Net weight in kg	0,20
Gross weight in kg	0,23