

E 975/2

Mobile unit

for the optimum cleaning of instruments on 2 levels.



- Equipped with an integrated spray arm
- Magnetic strip for automatic mobile unit recognition
- The vertical clearance for both levels is different

EAN: 4002514816449 / Material number: 07765790 / Old Material Number: 69597503D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	6
SPRI mesh trays [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	427
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	650
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	640
Net weight in kg	10,80
Gross weight in kg	12,80
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1