

APH 550 Sego  
Segolabel Miele Edition: Starter kit  
for producing sterile goods labels.



- Full package incl. printer, power unit, cables and 1000 labels
- Secure labelling of the packaged sterile goods
- Useful addition to Segosoft Miele Edition
- Useful addition for CUBE/CUBE X sterilisers
- Easy connection via USB cable

EAN: 4002514940502 / Material number: 07951400 / Old Material Number: 68814001D

<b>Product classification</b>	
Small sterilisers	•
<b>Machine affiliation</b>	
PST 1710	•
PST 1720	•
PST 2210	•
PST 2220	•
<b>Product category</b>	
Process documentation/operating documentation/communication	•
<b>Product properties</b>	
Colour	Grey
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	178
External dimensions, net width in mm	135
External dimensions, net depth in mm	240
External dimensions, gross height in mm	230
External dimensions, gross width in mm	300
External dimensions, gross depth in mm	240
Net weight in kg	2,45
Gross weight in kg	2,90
<b>Standards, certs, symbols</b>	
CE	•