

APH 596

Labels

for use with Segolabel printer APH 550 Sego (PRT 200).



- 1 roll of 1000 reusable adhesive labels
- Every label has second adhesive surface with a separating layer

EAN: 4002514940632 / Material number: 07951530 / Old Material Number: 68816201D

Product classification	
Small sterilisers	•
Large sterilisers	•
Machine affiliation	
PS 1201 B	•
PS 1201 B EXCELLENCE	•
PS 1202 B	•
PS 5662 V	•
PS 5666 V	•
PS 5669 V	•
PS 5962 H	•
PS 5966 H	•
PS 5968 H	•
PS 5969 H	•
Product category	
Process documentation/operating documentation/communication	•
Product properties	
Colour	White
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	55
External dimensions, gross width in mm	90
External dimensions, gross depth in mm	95
Net weight in kg	0,30
Gross weight in kg	0,39