

E 198

Insert

for optimum loading of up to 6 mesh and kidney trays.



- Suitable for 6 mesh or kidney trays
- 7 holders (6 compartments), height 160 mm mm, distance 67 mm mm

EAN: 4002514953373 / Material number: 07984040 / Old Material Number: 68519801D

Product category	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
Product category	
Insert for mesh trays	•
Insert for kidney dishes/basins	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	6
Distance between rails in mm	67
Rail height in mm	160
Area of application	
Reprocessing of dental instruments	•
Capacity	
Diverse mesh trays [number]	6
Kidney dishes [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	180
External dimensions, net depth in mm	495
External dimensions, gross height in mm	185
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	530
Net weight in kg	0,65
Gross weight in kg	0,90