


E 197

Insert

For optimum accommodation of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002514953380 / Material number: 07984090 / Old Material Number: 68519701D

Product classification	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	42
External dimensions, net width in mm	150
External dimensions, net depth in mm	225
External dimensions, gross height in mm	75
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	240
Net weight in kg	0,55
Gross weight in kg	0,62