

E 337/1

Insert

for the optimum loading of upright instruments.



- Suitable for probes, tweezers or mouth specula
- Equipped with 4 plastic holders
- 12 compartments approx. width 22 x 28 mm
- 4 compartments approx. width 25 x 28 mm
- 48 compartments approx. width 13 x 14 mm

EAN: 4002514953397 / Material number: 07984130 / Old Material Number: 69533702D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	148
External dimensions, net width in mm	173
External dimensions, net depth in mm	445
External dimensions, gross height in mm	205
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	490
Net weight in kg	1,91
Gross weight in kg	2,16