

E 802/1

Insert

for the optimum loading of different upright instruments.



- Suitable for loading probes, tweezers or mouth specula
- Equipped with 4 plastic holders
- 8 compartments approx. width 30 x 28 mm
- 8 compartments approx. width 16 x 28 mm
- 15 compartments approx. width 16 x 20 mm

EAN: 4002514953403 / Material number: 07984170 / Old Material Number: 69580202D

<b>Product classification</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8562	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
<b>Product category</b>	
Insert for instruments	•
Insert/module for dental instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	148
External dimensions, net width in mm	163
External dimensions, net depth in mm	295
External dimensions, gross height in mm	190
External dimensions, gross width in mm	195
External dimensions, gross depth in mm	430
Net weight in kg	0,65
Gross weight in kg	1,15