

E 806/1

Insert

for the optimum loading of 8 half-shell trays.



- 9 holders (8 compartments), W 295, spacing 37 mm mm
- Max. tray size 290 x 30 mm

EAN: 4002514969718 / Material number: 09013220 / Old Material Number: 69580603D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	8
Distance between rails in mm	37
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	114
External dimensions, net width in mm	305
External dimensions, net depth in mm	348
External dimensions, gross height in mm	125
External dimensions, gross width in mm	330
External dimensions, gross depth in mm	480
Net weight in kg	0,71
Gross weight in kg	0,84