

APH 591
Paper roll
for protocol printer APH 500 PRT100 and APH 510 PRT110.



- 5 paper rolls per pack
- Width: 76 mm mm

EAN: 4002515000038 / Material number: 09063410 / Old Material Number: 68817003D

| Product classification | |
|-------------------------------------------------------------|-------|
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Thermal disinfectors | • |
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Small sterilisers | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7825 | • |
| G 7826 | • |
| G 7831 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8581 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| PG 8593 | • |
| PS 1201 B | • |
| PS 1201 B EXCELLENCE | • |
| PS 1202 B | • |
| Product category | |
| Process documentation/operating documentation/communication | • |
| Product properties | |
| Colour | White |
| Area of application | |
| Reprocessing of pipettes/phials | • |
| Reprocessing of laboratory glassware | • |
| Reprocessing of slides | • |
| Reprocessing of butyrometers | • |
| Reprocessing of pharmacy/infusion bottles | • |
| Reprocessing of anaesthetic instruments | • |
| Reprocessing of MIS instruments | • |
| Reprocessing of ophthalmology instruments | • |
| Reprocessing of robotic instruments | • |
| Reprocessing of GYN instruments | • |
| Reprocessing of ENT instruments | • |
| Reprocessing of dental instruments | • |
| Reprocessing of theatre shoes | • |
| Reprocessing of basins | • |
| Reprocessing of kidney dishes | • |
| Reprocessing of baby bottles | • |
| Reprocessing of surgical instruments | • |
| Reprocessing of lumened instruments | • |

APH 591

Paper roll

for protocol printer APH 500 PRT100 and APH 510 PRT110.



EAN: 4002515000038 / Material number: 09063410 / Old Material Number: 68817003D

| | |
|---------------------------------------------------------|------|
| Reprocessing of orthopaedic surgery instruments | • |
| Reprocessing of intramedullary nails/bone marrow drills | • |
| Reprocessing of trays | • |
| Reprocessing of mesh trays | • |
| Reprocessing of measuring cylinders | • |
| Reprocessing of containers, e.g. soda lime containers | • |
| Reprocessing of different utensils | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 79 |
| External dimensions, net width in mm | 76 |
| External dimensions, net depth in mm | 79 |
| External dimensions, gross height in mm | 105 |
| External dimensions, gross width in mm | 165 |
| External dimensions, gross depth in mm | 210 |
| Net weight in kg | 1,39 |
| Gross weight in kg | 1,49 |