

KWT 6831 SG

Freestanding wine unit

with extra versatility and style, thanks to FlexiFrame and Miele Noteboard.



EAN: 4002515227084 / Material number: 09472250 / Old Material Number: 36683100AUS


User convenience	
Accessories required	XKS 3000 Z
Controls	
Independent temperature regulation of fridge and freezer zones	No
No. of temperature zones	3
Min. temperature setting range in °C	5
Max. temperature setting range in °C	20
Technical data	
Niche depth in mm	550
Appliance width in mm	700
Appliance height in mm	1920
Appliance depth in mm	746
Weight in kg	142,0
Wine conditioning zone in l	505
Total usable capacity in l	505
Storage time in case of fault in h	0
Voltage in V	220-240
Length of supply lead in m	2,6

- Simultaneous wine storage in separate temperature zones
- Odourless refrigeration section thanks to Active AirClean filters
- Present your wine in the perfect light – LED lighting
- Optimum wine protection through UV filters in tinted glass front
- Secure storage of wine bottles with pull-out racks


KWT 6831 SG

Freestanding wine unit


with extra versatility and style, thanks to FlexiFrame and Miele Noteboard.



KWT6831SG, KWT6832SGS, KWT6833SG,
KWT6834SGS (Installation drawings)



KWT6831SG, KWT6832SGS, KWT6833SG,
KWT6834SGS (Installation drawings)



KWT6831SG, KWT6833SG, KWT6832SGS,
KWT6834SGS (Installation drawings)