


U 503
Basket
For the optimum loading of cutlery.



- Close-meshed base
- Also suitable for sundae dishes
- Without cutlery inserts

EAN: 4002515250006 / Material number: 09508660 / Old Material Number: 67650301D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Basket for cutlery	•

U 503

Basket

For the optimum loading of cutlery.



EAN: 4002515250006 / Material number: 09508660 / Old Material Number: 67650301D

Basket for saucepans/bowls, etc.

Product properties

Material

Plastic

Colour

Green

Dimensions and weight

External dimensions, net height in mm

109

External dimensions, net width in mm

500

External dimensions, net depth in mm

500

External dimensions, gross height in mm

260

External dimensions, gross width in mm

520

External dimensions, gross depth in mm

530

Net weight in kg

1,90

Gross weight in kg

2,40