

HBB 71

Genuine Miele baking tray  
with PerfectClean finish.

EAN: 4002515261279 / Material number: 09519820 /  
Old Material Number: 22996249D



- Fits into FlexiClips or directly into the side runners
- Non-stick – no baking parchment or greasing required
- Especially easy cleaning thanks to Perfect Clean finish
- For Generation H 2000/6000/7000, DGC 6000/7000, DGC HC Pro 7000
- Dimensions:  
44.8 cm x 38.6 cm x 2.2 cm