

HUBB 91
Genuine Miele multi-purpose tray
with PerfectClean finish.

EAN: 4002515261354 / Material number: 09520680 / Old Material Number: 22996257D



- As a drip tray for meat juices/fat from roasting and grilling
- For use when roasting with the grilling and roasting insert
- Fits into FlexiClips or directly into side racks
- For 90 cm wide Generation H 2xxx/H 6xxx/H 7xxx appliances
- Dimensions: 72 cm x 36.55 cm x 5.6 cm
- Especially easy cleaning thanks to Perfect Clean finish