

APH 332
 Adapter cable
 for the connection of W&H Lisa 317, 322, 517, 522.

EAN: 4002515350959 / Material number: 09553890 / Old Material Number: 68812503D

Product category	
Small sterilisers	•
Product category	
Process documentation/operating documentation/ communication	•
Area of application	
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of mesh trays	•
Reprocessing of containers, e.g. soda lime containers	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0

- In addition to connection cable type 2
- Length 3 m
- For connecting the machine to Segosoft Miele Edition