


A 303
Module
for the optimum loading of up to 98 pipettes.



- Suitable for 98 measuring and full pipettes, for example
- Height of holding frame 150 mm
- 1 unit can be inserted into each of the baskets A 150 and APLW 106
- Vertical clearance: max. 450 mm (in A 150) and 605 mm (in APLW 106)

EAN: 4002515460825 / Material number: 09862560 /
Old Material Number: 69630300D

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Machine affiliation	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
Product category	
Insert/module for pipettes	•
Area of application	
Reprocessing of pipettes/phials	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
Pipettes [number]	98
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	185
External dimensions, net width in mm	222
External dimensions, net depth in mm	471
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	2,50
Gross weight in kg	3,50