

A 800

Tubular filter

For the safe and economic reprocessing of hollow instruments.



- Filter in ply mesh stainless steel mesh for mobile units
- Mesh size 63 µm
- Reprocessable with cleaning brush A 804

EAN: 4002515538081 / Material number: 09960680 / Old Material Number: 68680000D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Filter	•
Product properties	
Material	Stainless steel
Dimensions and weight	
External dimensions, net height in mm	23
External dimensions, net width in mm	23
External dimensions, net depth in mm	405
External dimensions, gross height in mm	45
External dimensions, gross width in mm	220
External dimensions, gross depth in mm	445
Diameter in mm	23
Net weight in kg	0,36
Gross weight in kg	0,88