

PG 8596
 Dispensing cabinet
 For process chemicals and dispenser modules



- 2 removable drawers for 3 x 5 l canisters each
- Internal dimensions H 285-690 , W 250, D 425-555 mm
- Each level 2 holding options per level for dispenser modules
- Drip tray facilitates better residual emptying of the canister

EAN: 4002515558324 / Material number: 10087340 /
 Old Material Number: 69859602D

Product category	
Laboratory washers	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Dispenser system	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	835
External dimensions, net width in mm	300
External dimensions, net depth in mm	600
External dimensions, gross height in mm	900
External dimensions, gross width in mm	400
External dimensions, gross depth in mm	740
Net weight in kg	26,50
Gross weight in kg	30,00