

## UBS 2

### Conversion kit

for use with selected series G 78 mobile units in PG 85 machines.



- For older injector units without drying unit connection

EAN: 4002515588987 / Material number: 10131100 / Old Material Number: 69215015D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8593	•
<b>Product category</b>	
Conversion kit	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	90
External dimensions, net width in mm	129
External dimensions, net depth in mm	295
External dimensions, gross height in mm	125
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	350
Net weight in kg	0,51
Gross weight in kg	0,71