

UBS 3

Conversion kit

For use with selected series G 78 upper baskets in PG 85 machines.

EAN: 4002515588994 / Material number: 10131110 / Old Material Number: 69215016D



- For older load carriers without dryer connection
- For older load carriers with welded water supply

| | |
|---|------------------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Thermal disinfectors | • |
| Machine affiliation | |
| PG 8581 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| PG 8593 | • |
| Product category | |
| Conversion kit | • |
| Product properties | |
| Colour | Stainless Steel, Black |
| Dimensions and weight | |
| External dimensions, net height in mm | 90 |
| External dimensions, net width in mm | 305 |
| External dimensions, net depth in mm | 176 |
| External dimensions, gross height in mm | 100 |
| External dimensions, gross width in mm | 310 |
| External dimensions, gross depth in mm | 210 |
| Net weight in kg | 0,59 |
| Gross weight in kg | 0,69 |