

A 602 Module

For the optimum reprocessing of instruments used in robot-assisted surgery.



- Suitable for 8 instruments, latest generation (Xi)

EAN: 4002515585030 / Material number: 10131930 / Old Material Number: 69660201D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Module/rack for robotic instruments	•
Area of application	
Reprocessing of robotic instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Robotics instruments [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	303
External dimensions, net width in mm	555
External dimensions, net depth in mm	318
External dimensions, gross height in mm	540
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	500
Net weight in kg	6,30
Gross weight in kg	7,80