


A 807
Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments, surgical clamps, older generation, S/Si

EAN: 4002515585078 / Material number: 10131980 / Old Material Number: 69680701D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Module/rack for robotic instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of robotic instruments	•
Capacity	
Robotics instruments [number]	1
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	77
External dimensions, net width in mm	41
External dimensions, net depth in mm	165
External dimensions, gross height in mm	25
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	485
Net weight in kg	0,37
Gross weight in kg	0,49