

MC/2

Transfer trolley

height adjustable, for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- 2 loading racks (trough shaped) as well as handle and cover plate

EAN: 4002515588628 / Material number: 10138690 / Old Material Number: 69200109D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8535 | • |
| PG 8536 | • |
| PG 8581 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| PG 8593 | • |
| Product category | |
| Storage/transport system | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Dimensions and weight | |
| External dimensions, net height in mm | 1003 |
| External dimensions, net width in mm | 600 |
| External dimensions, net depth in mm | 805 |
| External dimensions, gross height in mm | 1119 |
| External dimensions, gross width in mm | 835 |
| External dimensions, gross depth in mm | 1280 |
| Net weight in kg | 34,20 |
| Gross weight in kg | 52,40 |