


KWT 6833 SG

Freestanding wine unit

with extra versatility and style, thanks to FlexiFrame and Miele Noteboard.



- Simultaneous wine storage in separate temperature zones
- Odourless refrigeration section thanks to Active AirClean filters
- Present your wine in the perfect light – LED lighting
- Optimum wine protection through UV filters in tinted glass front
- Secure storage of wine bottles with pull-out racks

EAN: 4002515716663 / Material number: 10409080 / Old Material Number: 36683300EU1

Appliance category	
Wine conditioning unit	•
Construction type	
Freestanding appliance	•
Door hinge	right
Convertible door hinging	•
Design	
Casing colour	Stainless Steel
Wine compartment lighting	LED illumination
User convenience	
Networking with Miele@home	•
Accessories required	XKM 3100 W
DynaCool	•
SoftClose	•
Concealed door opener	Recessed grip
No. of wooden racks	13
Number of FlexiFrame racks	10
NoteBoard	Yes
Bottle presenter	•
Control system	
Operating concept	TouchControl
Fridge compartment can be switched off	No
Independent temperature regulation of fridge and freezer zones	No
No. of temperature zones	3
Min. temperature setting range in °C	5
Max. temperature setting range in °C	20
Efficiency and sustainability	
Energy efficiency class (A–G)	G
Annual energy consumption in kWh	205,50
Energy consumption in 24 h in kWh	0,56
Safety	
Lock function	•
Acoustic door alarm	•
Acoustic temperature alarm	•
Optical door alarm	•
Optical temperature alarm	•
Technical data	
Appliance width in mm	700
Appliance height in mm	1920
Appliance depth in mm	746
Weight in kg	142,0
Climate class	SN-ST
Refrigerator compartment in l	503
Wine compartment in l	503
Total usable capacity in l	503
0.75 l Bordeaux bottle capacity (no. of bottles)	178 bottles
Storage time in case of fault in h	0
Acoustic airborne noise emission class (A–D)	B
Sound emissions in dB(A) re 1 pW	35

KWT 6833 SG

Freestanding wine unit

with extra versatility and style, thanks to FlexiFrame and Miele Noteboard.



EAN: 4002515716663 / Material number: 10409080 / Old Material Number: 36683300EU1

Current consumption (mA)	1500
Voltage in V	220-240
Fuse rating in A	10
Number of phases	1
Frequency in Hz	50
Length of supply lead in m	2,8

Accessories included


Active AirClean filter

•


KWT 6833 SG

Freestanding wine unit


with extra versatility and style, thanks to FlexiFrame and Miele Noteboard.



KWT6831SG, KWT6832SGS, KWT6833SG,
KWT6834SGS (Installation drawings)



KWT6831SG, KWT6832SGS, KWT6833SG,
KWT6834SGS (Installation drawings)



KWT6831SG, KWT6833SG, KWT6832SGS,
KWT6834SGS (Installation drawings)