

A 310  
Insert  
for the optimum loading of 8 theatre shoes.



- For use in upper baskets
- For loading 8 theatre shoes up to size 41
- Can be used in upper baskets A 101 and A 102

EAN: 4002515816349 / Material number: 10632610 / Old Material Number: 69631000D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for theatre shoes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Theatre shoes [number]	8
<b>Dimensions and weight</b>	
External dimensions, net height in mm	75
External dimensions, net width in mm	405
External dimensions, net depth in mm	435
External dimensions, gross height in mm	100
External dimensions, gross width in mm	425
External dimensions, gross depth in mm	450
Net weight in kg	0,50
Gross weight in kg	0,83