

A 817

Irrigation sleeves

for optimum cleaning of ENT instruments, for example.



- Suitable for mobile injector units and baskets, screw connection
- Including cap and mounting clip

EAN: 4002515880593 / Material number: 10706260 / Old Material Number: 69681700D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8534	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Diameter in mm	24
Net weight in kg	0,02
Gross weight in kg	0,02