

A 819

Irrigation sleeves

for optimum cleaning of lumened instruments.

EAN: 4002515990278 / Material number: 10721750 / Old Material Number: 68681900D



- E.g. for magnetostrictive ultrasound scaler tips, e.g. Cavitron
- Int. dimension 108 mm
- Including cap with an opening of Ø 10 mm
- With mounting clip

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
Product category	
Holder	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of dental instruments	•
Capacity	
Ultrasonic tips [number]	1
Dimensions and weight	
External dimensions, net height in mm	10
External dimensions, net width in mm	100
External dimensions, net depth in mm	114
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0,02
Gross weight in kg	0,02