

TC 900

Dresser panelling

For machines 900 mm wide without a steam condenser or drying unit.



- Includes service door with lock
- Diaphragm wall panelling up to ceiling must be provided on site

EAN: 4002515991756 / Material number: 10880700 / Old Material Number: 69252003D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
<b>Product category</b>	
Cover	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	840
External dimensions, net width in mm	845
External dimensions, net depth in mm	1045
External dimensions, gross height in mm	990
External dimensions, gross width in mm	845
External dimensions, gross depth in mm	1200
Net weight in kg	48,90
Gross weight in kg	63,90