

DS 10

Membrane pump

For 10-litre containers, for dispensing liquid media.



- For retrofitting additional dispenser pumps
- For dispensing detergents, rinsing agents and neutralising agents
- Incl. siphon for 10-litre media container
- Scope of delivery: Dispenser pump, lance, volume control & hoses

EAN: 4002515992241 / Material number: 10880870 / Old Material Number: 69249101D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Dispenser system	•
Standard electrical connection	
Number of phases	0
AC/DC	DC
Voltage in V	24
Frequency in Hz	0
Total rated load in kW	0,01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	200
External dimensions, gross width in mm	480
External dimensions, gross depth in mm	150
Net weight in kg	1,05
Gross weight in kg	1,25