


RP 900

Plinth, on wheels

For machines 900 mm wide, incl. drip tray without openings for floor drain.



- Suitable for single-door machines, 900 mm wide
- Castors allow the machine to be pulled forward for servicing
- Only suitable for machines with a drain pump
- Easy access for servicing

EAN: 4002516032465 / Material number: 10941940 / Old Material Number: 69253007D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
Product category	
Plinth/base	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	100
External dimensions, net width in mm	900
External dimensions, net depth in mm	990
External dimensions, gross height in mm	155
External dimensions, gross width in mm	940
External dimensions, gross depth in mm	1200
Net weight in kg	30,80
Gross weight in kg	32,00