

A 503

Mobile unit

no levels, for cleaning laboratory glassware/utensils.



- For holding inserts and mesh trays
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516050445 / Material number: 10973900 / Old Material Number: 69650301D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit with module connection	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	433
External dimensions, net width in mm	790
External dimensions, net depth in mm	650
External dimensions, gross height in mm	790
External dimensions, gross width in mm	860
External dimensions, gross depth in mm	715
Net weight in kg	17,00
Gross weight in kg	20,00