

A 314

Mesh trays

For the optimum reprocessing of ophthalmic lumened instruments.



- Incl. 1 injector manifold to connect ophthalmic instruments
- Mesh tray with removable lid
- Particularly suitable for A 207 mobile unit

EAN: 4002516321262 / Material number: 10979670 / Old Material Number: 69631400D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of ophthalmology instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	38
External dimensions, net width in mm	223
External dimensions, net depth in mm	151
External dimensions, gross height in mm	75
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	245
Net weight in kg	0,72
Gross weight in kg	0,85