


A 613
Module
for the optimum reprocessing of volumetric and measuring pipettes.



- Equipped with 98 injector nozzles
- Suitable for measuring pipettes with a volume of 1–50 ml
- Suitable for volumetric pipettes with a volume of 1–100 ml
- For use in A 503 mobile unit

EAN: 4002516063360 / Material number: 10999070 / Old Material Number: 69661300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Capacity	
Pipettes [number]	98
Dimensions and weight	
External dimensions, net height in mm	194
External dimensions, net width in mm	254
External dimensions, net depth in mm	604
External dimensions, gross height in mm	200
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	610
Net weight in kg	4,40
Gross weight in kg	6,40