


A 623  
Module  
for the optimum loading of measuring cylinders.



- Equipped with 4 nozzles
- For measuring cylinders with 1–2 l capacity
- Pivoting supports for ease of use

EAN: 4002516071440 / Material number: 11025560 / Old Material Number: 69662300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for measuring cylinder	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Capacity	
Measuring cylinders [number]	4
Dimensions and weight	
External dimensions, net height in mm	437
External dimensions, net width in mm	255
External dimensions, net depth in mm	608
External dimensions, gross height in mm	540
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	500
Net weight in kg	6,30
Gross weight in kg	7,80