

APH 592
Colour transfer belt
for protocol printer APH 510 PRT110.



EAN: 8809166675907 / Material number: 11053500 / Old Material Number: 68817005D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PWD 8541	•
PWD 8545	•
PWD 8546	•
PWD 8549	•
Product category	
Process documentation/operating documentation/ communication	•
Product properties	
Colour	Red, Black
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of slides	•
Reprocessing of butyrometers	•
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of measuring cylinders	•
Reprocessing of containers, e.g. soda lime containers	•
Reprocessing of bedpans and urine bottles	•

APH 592
Colour transfer belt
for protocol printer APH 510 PRT110.



EAN: 8809166675907 / Material number: 11053500 / Old Material Number: 68817005D

Reprocessing of different utensils

•

Dimensions and weight

External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	55
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	250
Net weight in kg	0,04
Gross weight in kg	0,05

Standards, test certificates and symbols

CE

•