


A 837

Irrigation sleeve set

For reprocessing lumened instruments used in ENT applications in particular.



- For use in A 105/1 upper basket
- Incl. cap and mounting clip
- Pack of 22 irrigation sleeves (A 817)

EAN: 4002516204107 / Material number: 11054690 / Old Material Number: 69783700D

Product classification	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Irrigation sleeve	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of ENT instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	55
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	250
Net weight in kg	0,53
Gross weight in kg	0,53