


A 302/2

Module

for optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 32 injector nozzles
- Part of the EasyLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth
- Particularly suitable for A 844 and A 845 (nozzle diameter 2.5 mm)

EAN: 4002516087656 / Material number: 11056790 / Old Material Number: 69630202D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	32
Standard electrical connection	
Number of phases	0
Dimensions and weight	
External dimensions, net height in mm	72
External dimensions, net width in mm	195
External dimensions, net depth in mm	479
External dimensions, gross height in mm	200
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	1,09
Gross weight in kg	2,16