

A 861  
Holder rack  
for A 301/4 and A 301/5 injector modules.



- Arrangement: 3 x 6 fields
- For use on A 301/4 and A 301/5 injector modules
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516087670 / Material number: 11056810 / Old Material Number: 69686100D

<b>Product category</b>	
Laboratory washers	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	18
Compartment width in mm	74
Compartment length in mm	74
Part of the EasyLoad system	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	80
External dimensions, net width in mm	200
External dimensions, net depth in mm	350
External dimensions, gross height in mm	100
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	495
Net weight in kg	0,10
Gross weight in kg	1,88