

A 861
Holder rack
for injector modules A 301/4 and A 301/5



- Arrangement: 3 x 6 fields
- For use on A 301/4 and A 301/5 injector modules
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516087670 / Material number: 11056810 / Old Material Number: 69686100D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	18
Compartment width in mm	74
Compartment length in mm	74
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	80
External dimensions, net width in mm	200
External dimensions, net depth in mm	350
External dimensions, gross height in mm	100
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	495
Net weight in kg	0,10
Gross weight in kg	1,87