

APLW 052/1

Set of injector nozzles

With plastic supports, Ø 4 mm, length 175 mm.

EAN: 4002516130918 / Material number: 11094000 / Old Material Number: 69305201D



- Suitable for e.g. lower basket APLW 095
- Packaging unit with 39 injector nozzles
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 5 mm

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Dimensions and weight	
External dimensions, gross height in mm	190
External dimensions, gross width in mm	310
External dimensions, gross depth in mm	310
Net weight in kg	0,78
Gross weight in kg	1,05