

APLW 001

Mobile unit

for optimum loading of laboratory glassware on the lowest level.



- Maximum loading width: 470 mm
- Maximum loading depth: 490 mm

EAN: 4002516115298 / Material number: 11104480 / Old Material Number: 69300100D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for beakers	•
Mobile unit for inserts	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	135
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	190
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	4,20
Gross weight in kg	5,90